



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Respirator PURITAN BENNETT 840 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Made in U.S.A.,
Power supply: 230 V,
Emergency power: via 2 new built-in batteries 12 V / 7 Ah,
Frequency: 50/60 Hz,
Ventilator system
Puritan Bennett 840
is equipped with an intuitive touch screen that facilitates assessment of patient data and control parameters,
The SmartAlert system increases safety by capturing alarms and providing useful information. The system supports various modes and types of ventilation, providing flexible patient care,
Features:
Touch screen enables intuitive control and monitoring,
SmartAlert system for proactive alarm management,
Flexible modes and types of breath delivery adapted to patient needs,
Support for different patient groups (children, adults),
Built-in PURITAN 806 compressor,
Specification:
Ventilation modes: A/C, SIMV, SPONT, BILEVEL,
Ventilation type: Invasive / Non-invasive,
Breath forced: PC / VC / VC+,
Trigger mode: V-TRIG / P-TRIG,
Pressure support (PSUPP): 0 to 70 cm H₂O,
Expiratory sensitivity range (ESENS): 1% to 80%;
Ventilation rate of tidal volume without apnea (VT): 25 to 2500 ml,
Breathing frequency without apnea (f): 2.0 to 40 / min,
Peak inspiratory flow (V.MAX):
3.0 to 60 l / min, for pediatric patient,
3.0 to 150 l / min, for adult patient,
Inspiratory pressure (PI): 5 to 90 cm H₂O,
Inspiration time without apnea (TI): 0.2 to 8.0 seconds,
Apnea interval (TA): 10 to 60 seconds,
Apnea expiration time (TE): TE ≥ 0.2 seconds,
End inspiratory pressure: PEEP: - 20 to 130 cm H₂O,
Built-in gas connections: AIR, O₂,
Menu language: Polish,
Ventilator placed on a dedicated mobile cart,
Included in the set:
User manual in English (PDF),
Power cord,
2x new bacterial filter,
New bacterial exhaust filter,
Overall dimensions: 60 x 65 x 148 cm,
Weight: 100 KG,
Has current inspection and is ready for work,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months on domestic market (Poland),
3 months on international market,
Option to extend warranty for additional charge up to 12 months on domestic and international markets,
Financing option (Poland only): Installments, Leasing, Loan,

