



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MAQUET MEDAP TWISTA SILENT POWER 1070 surgical vacuum (Reconditioned)**

Refurbished (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
German production, (2017),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
MEDAP TWISTA SP 1070  
to exceptionally efficient surgical aspirator designed specifically for reliable, continuous operation in the operating room,  
Technical specification:  
Suction performance:  
50 Hz:  $58 \pm 6$  l/min.,  
60 Hz:  $68 \pm 7$  l/min.,  
Max. vacuum  
90 kPa,  
900 mbar,  
Vacuum regulation range: From 0 to -1000 mbar,  
Operating mode: Continuous operation,  
Noise level: 53 dB(A),  
Included:  
Power cord,  
New autoclavable Serres 2 L jar with a new disposable insert,  
Dimensions: 50 x 60 x 90 cm,  
Weight: 26 KG,  
Has a current inspection and is ready to operate,  
Technical Passport issued (Service Report) valid for 12 months,  
Warranty:  
6 months on the domestic market (Poland),  
3 months on the international market,  
Warranty extension available for an additional fee up to 12 months for both domestic and international markets,  
Finance options (Poland only): Installments, Leasing, Loan,

