



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **PROGETTI RESCUE LIFE defibrillator (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Italian production,
- Power supply: 230 V,
- Emergency power supply: via new built-in NiMh battery,
- Frequency: 50/60 Hz,
- **The PROGETTI RESCUE LIFE** is a professional manual clinical defibrillator designed for advanced resuscitation procedures performed by trained personnel. Standard equipment includes an SpO2 module and AED functionality.
- **Technical specifications :**
  - Screen: large 7" color TFT display,
  - ECG signal: from paddles or 10-lead cable,
  - Heart rate,
  - **ECG monitor:**
    - Travel: 5 / 12.5 / 25 / 50 mm/sec,
    - Amplitude (feature) 5 / 10 / 20 mm/mV,
    - Heart rate – 20 to 300 bpm ( $\pm 2\%$ ),
    - Alarm: HR alarm setting in the range of 20-250 bpm,
  - **Pulse type:**
    - BTE (biphasic exponential pulse); compensated impedance,
    - Energy: Manual: 1-230J,
    - Energy: AED: 150J – constant,
    - HiCAP (Large Storage Capacitor) technology,
    - Defibrillation energy selection: 1 – 230J – multifunctional knob on the panel,
  - **Charging time:**
    - 100 discharges with charged batteries,
    - Shock synchronization mode with R (cardioversion),
    - Charging/defibrillation with a button on the paddles,
  - **Memory: Internal 1 Gb flash memory, up to 300 hours of ECG recording, events, review,**
  - **Printer:**
    - Integrated thermal printer with a paper width of 58 mm for ECG waveforms and events, printout documentation with values,
  - **Power supply:**
    - High-capacity NiMh battery with short charging time (95% in 2 hours),
    - Mains power/charging: 100-230 V,  $\pm 10\%$  50/60 Hz,
- **Menu in: English,**
- **The device has visible signs of wear that do not affect the comfort of use,**
- **Included in the set:**
  - PROGETTI pediatric spoons,
  - New 10-lead ECG cable,
  - User manual in English (PDF),
  - Power cord,
- Dimensions: 37 x 34 x 34 cm,
- Weight: 5.5 KG,
- Cat D154
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months to the international market,**
  - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

