



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **VERMEIREN CLUB hospital bed (Reconditioned)**

Refurbished (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Belgian production, (2025),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
New mattress with medical certificate in navy blue,  
VERMEIREN VARIO CLUB hospital bed  
is a basic care bed intended for hospital wards, rehabilitation and also home care,  
High quality and functional standard make  
VARIO CLUB  
an excellent bed supporting the care and rehabilitation process and improving the quality of care for lying patients,  
Specification:  
Four-segment lying surface,  
Smooth adjustable bed height (40 to 80 cm) via remote,  
Anti-Trendelenburg position,  
Smooth angle adjustment of the headboard via remote,  
Smooth angle adjustment of the thigh segment via remote,  
LINAK actuators,  
Beige beech color housing,  
Overhead lift and rails included,  
In the set:  
User manual in FR, ND, GER, IT, ESP, PL, SK,  
Patient overhead lift,  
Dimensions:  
Mattress: 90 x 196 cm,  
Overall: 103 x 218 x 84/124 cm,  
Maximum load: 165 KG,  
Weight: 62 KG  
Has current inspection and is ready for operation,  
Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months on the domestic market (Poland),  
3 months on the international market,  
Possibility to extend warranty for an additional fee up to 12 months for domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan,

