



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Stretcher trolley for transporting patients MERIVAARA EMERGO (Reconditioned)

Refurbished (used),

Technical condition: very good,

Visual condition: good,

Real product photos,

Finnish production

Medical certificate mattress in black,

Stretcher trolley

MERIVAARA EMERGO

is a device for transport and care, designed with hospital needs in mind. Its modular construction allows operation according to clinical needs, whether for patient transport, emergencies, intensive care units, convalescence phase, or even during minor procedures,

The cot frame combines sturdiness, safety and comfort, guaranteeing reliable operation even in the most demanding hospital conditions, EMERGO

stands out with exceptional mobility. Large, anti-static wheels with a diameter of 150 mm with a central lock provide smooth and safe movement on all kinds of hospital floors,

Specification:

Manual adjustment:

Height adjustable using foot manipulator in the range: 63 - 100 cm,

Backrest tilt range: 70°,

Trendelenburg tilt: 14°

Anti-Trendelenburg tilt: 7°,

Double retractable rails,

Built-in bumpers,

Built-in central brake,

Dimensions:

Mattress: 65 x 192 cm,

Overall: 78 x 215 x 63/100 cm,

Maximum load: 225 KG,

Weight: 85 KG,

Has current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for domestic market (Poland),

3 months for international market,

Possibility to extend warranty for a fee up to 12 months for domestic and international markets,

Financing available (Poland only): Installments, Leasing, Loan,

