



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Polometer (perimeter) OPTOPOL PTS-1000 (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Polish production, (2011),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Software: 4.4.3,
- **The PTS-1000 perimeter** is a device designed for visual field testing. It is an essential piece of equipment in every ophthalmologist's office. It enables quick and precise testing using a static stimulus,
- Offers static and kinetic stimuli with all Goldmann stimulus sizes and all stimulus colors used in perimetry,
- A wide range of analytical strategies and tools combined with an intuitive interface provide a powerful diagnostic tool,
- Built-in maps and the ability to create your own allow you to explore only the areas you need,
- **Specification:**
 - Perimetry type: Static / kinetic,
 - Ball radius: 300 mm,
 - Field range: 100° (up to 160° with horizontal fixation shift),
 - Goldmann stimulus size: I - V,
 - Stimulus colors: white / blue / red / green,
 - Fixation method: blind spot position control or digital pupil detection,
 - Maximum stimulation intensity,
 - Stimulus velocity in kinetic perimetry: 1000 asb,
 - Patient response time: From 0.1 to 9.9 seconds,
 - **Strategies:**
 - Kinetic,
 - Screen,
 - Threshold,
 - Fast threshold,
 - 3-zone,
 - BSV,
 - FLICKER,
 - SWAP,
 - BDT,
- **Equipment:**
 - **AIO Computer ACL OR-PC 19 LP:**
 - **ACL** is a medical, fanless, compact PC,
 - Designed for continuous 24/7 operation in the patient environment,
 - Suitable for convection cooling,
 - Equipped with an internal 19" TFT display with a wide viewing angle,
 - It has an easy-to-clean construction,
 - **Screen:**
 - Resolution: 1280 x 1024 px,
 - Contrast 1500:1,
 - Brightness 300 cd/m²,
 - Viewing angle: 178°/178°,
 - Operating system: Windows 7 Professional,
 - RAM: 4 GB,
 - Processor: i5-3337U (1.80 GHz),
 - Disk: 128 GB.
 - **It has the following inputs and outputs:**
 - HDMI,
 - DVI,
 - LAN – 2 pcs,
 - USB 2.0 – 4 pcs,
 - COM 1 and 2,
 - SOUND,
- **The equipment is placed on a mobile table MDT MD-2N with electric height adjustment in the range of 76 - 96 cm,**
- **Included in the set:**
 - Power cord,
 - Patient button,
 - DELL Wired Mouse,
 - DELL Wired Keyboard,
 - LEXMARK MS312dn printer,
 - Patient eye cover,
- Overall dimensions: 100 x 50 x 148 cm,
- Weight: 90 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
 - **The warranty period does not cover light bulbs and other light sources installed in the device,**
- Financing options (Poland only): Installments, Leasing, Loan,

