



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Dental unit suction external GREELOY GS-E1000 (Reconditioned)

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 1000 W,  
UV lamp: 6W,  
Intelligent CNC panel, fully touch controls, you can adjust the device output with just one click; the digital knob displays operating parameters, it is easy and convenient,  
Specification:  
Control mode: touchscreen - two adjustable speeds,  
Suction flow: 4000 l/min.,  
Filtration efficiency: 99.99%,  
Noise:  $\leq 60$  dB,  
Primary filtration: 5-25 $\mu$ m,  
Secondary filtration: 0.3 $\mu$ m,  
Free adjustment of the working arm within lateral rotation and height,  
In the set:  
Power cable,  
Dimensions: 34 x 34 x 180 cm,  
Weight: 30 KG,  
Device equipped with new filters,  
Has current inspection and is ready for use,  
Technical Passport issued (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Warranty extension available for an additional charge up to 12 months for both domestic and international markets,  
Warranty period does not cover bulbs and other light sources installed in the device,  
Financing available (Poland only): Installments, Leasing, Loan,

