



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STRYKER 721 Medical Transport Trolley (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: good,
Actual product photos,
Made in USA,
New mattress with medical certification in black,
STRYKER 721
is a transport bed with a simple construction that provides safe and efficient transport,
Used mainly for general purposes, transport, procedures and convalescence. Also performs well in clinics and other healthcare facilities,
Manual adjustment:
Using foot switches:
Bed height range: 64 - 98 cm,
Trendelenburg position: +18°,
Anti-Trendelenburg position: -18°,
Backrest adjustment with gas spring range: 0° - 72°,
Equipped with:
Folding side rails,
Corner roll-off bumpers,
Central brake and steering system that provide mobility for easier traction and obstacle avoidance.
Four large wheels with a diameter of 20 cm,
IV stand height 122 cm,
Dimensions:
Mattress: 58 x 188 cm,
Overall: 70 x 214 x 64 - 98 cm,
Maximum load: 228 KG,
Weight: 100 KG,
Has a current inspection and ready for work,
Technical Passport (Service Report) issued and valid for 12 months,
Warranty:
6 months on domestic market (Poland),
3 months on international market,
Warranty extension available for an extra charge up to 12 months for domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

