



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



MEDTRONIC LIFEPAK 20E defibrillator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real product photos,
- Made in the USA (2015)
- Power supply: 230 V,
- Emergency power supply: via 1 new built-in battery,
- Frequency: 50/60 Hz,
- Power: 120 VA,
- The Lifepak 20E is a biphasic, portable defibrillator that provides a complete response system for acute cardiac emergencies.
- Technical specifications:
 - AED mode,
 - Synchronization,
 - Monitoring ECG parameters,
 - Easy-to-read 6.2" display
- Built-in printer (paper width: 50 mm),
- Possibility to select menu in: ENG, GER, FR, ESP, SWE, IT, DUN, FI, NED, NOR, PL, PT, BR, HUN,
- Included in the set:
 - User manual in Polish (PDF),
 - 3-lead ECG cable,
 - Power cord,
 - QUICK-COMBO cable,
- Dimensions: 26 x 24 x 21 cm,
- Weight: 6.5 KG,
- CAT D155,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months to the international market,
 - Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,
- Financing options (Poland only): Installments, Leasing, Loan,

