



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HD-LED treatment lamp BRANDON ASTRALITE ALE10CJ (ceiling) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- UK Production, (2017),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- **The Astralite lamp** is designed so that the lamp head can be positioned around the table to examine all parts of the patient's body. The freely adjustable lamp head allows for easy adjustment of the light beam to any angle,
- **Characteristics:**
 - Antibacterial additives destroy harmful bacteria. Smooth, continuous surfaces are easy to clean, and parting lines are minimized. The Astralite lamp head is sealed to prevent the ingress of dirt and fluids, preventing contamination.
 - Biomaster antibacterial additive is injected into the plastic parts of the lamp head, and "Ultra Fresh" additive is incorporated into the paint on the lamp head. Both additives have antibacterial properties,
 - The light color can be adjusted to suit the skin, tissue or wound type to select optimal illumination,
 - The lamp produces a large light field with a maximum diameter of 250 mm and a working distance of 70-140 cm.
- **Specification:**
 - Maximum light intensity - 100,000 lux.,
 - Power consumption -
 - Color temperature - 3500 - 5000°K,
 - Field size - 250 mm.,
 - Beam size at 1 m.- 180 mm.,
 - Color rendering index Ra: ≥ 95 ,
 - Color rendering (red) R9: ≥ 95 ,
 - Lighting depth - 600 mm.,
- **The minimum ceiling height for installing the lamp is: 300 / 310 cm,**
- **Included in the set:**
 - Column for hanging a lamp,
 - Cover,
 - Mounting kit,
- Dimensions:
 - Canopy: 30 cm,
 - 2-segment arm spacing: 156 cm,
 - Mounting column: 60 cm,
- Weight: 24 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
 - **The warranty period does not cover light bulbs and other light sources installed in the device ,**
- Financing options (Poland only): Installments, Leasing, Loan,

