



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



### **Patient monitor/kardiometer MINDRAY PM-9000 (Reconditioned)**

Refurbished (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
German production,  
Power supply: 230 V,  
Backup power: via 1 new built-in 12 V / 2 Ah battery,  
Frequency: 50/60 Hz,  
Power: 140 VA,  
MINDRAY PM-9000  
is a cardiac monitor monitoring patient vital functions,  
Specification:  
Color TFT display with 12.1" diagonal,  
Control by navigation knob and built-in buttons,  
Maximum 96-hour graphical and tabular trends of all parameters,  
Built-in measurement connectors:  
SPO2,  
ECG,  
NIBP,  
Temperature measurement [T1] | [T2],  
Menu in Polish,  
In the set:  
Power cable,  
New NIBP cuff with air hose,  
New SPO2 clip sensor,  
ECG Lead Cable 3-lead,  
Dimensions: 31 x 13 x 27 cm,  
Weight: 4.7 KG,  
Has current service and is ready for use,  
Technical Passport issued (Service Report) valid for 12 months,  
Warranty:  
6 months on domestic market (Poland),  
3 months on international market,  
Option to extend warranty for a surcharge to 12 months on domestic and international markets,  
Financing option (Poland only): Installments, Leasing, Loan,

