



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS KV-5 endoscopic suction machine (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Produced in the U.K.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 130 VA,
Simple and intuitive operation,
Suitable for all endoscopic procedures,
Variable suction with vacuum gauge: from 0 to -1 bar,
Maximum suction pressure: -85 kPa,
Quiet pump operation,
In the set:
Power cable,
User manual in English (PDF),
New Seeres jar with disposable insert of 1000 ml capacity,
New hydrophobic microbiological filter,
Dimensions: 27 x 39 x 22 cm,
Weight: 13 KG,
Cat 02,
Has a current inspection and is ready to operate,
Technical Passport issued (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee up to 12 months for both domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

