



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Video processor- Olympus CV-170 light source (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Manufactured in Japan,
Power: 230 V,
Frequency: 50/60 Hz,
Power: 200 VA,
OLYMPUS CV-170
is an integrated video processor with LED light source, intended for high-quality endoscopic diagnostics in medical offices and small treatment rooms,
Offers image in HDTV resolution
HDTV
and technology
NBI (Narrow Band Imaging)
, which support early detection of cancerous changes, especially in otolaryngology (laryngology, rhinoscopy) and urology (cystoscopy),
Technical specification:
HDTV image quality: Provides a clear full-screen image, allowing observation of mucosal structures and capillary vessels,
NBI technology: Improves visibility of surface vascularization, which is crucial in early diagnosis of precancerous and cancerous lesions,
LED light source,
Image adjustments:
White balance (front panel button),
Color tint adjustment (red, blue, chroma) in ± 8 degrees,
Autodyfragn mode selection (Peak or Average),
Image enhancement settings (structure enhancement or edge enhancement) at 3 levels,
Automatic gain control (AGC) and noise reduction,
Pre-freeze function: Automatically selects the sharpest frame from the sequence before pressing the "freeze" button,
Records up to 50 patient data sets in advance,
Ability to save up to 20 user profiles,
Menu language: Polish,
Ability to equip with Olympus MAJ-2004 / MAJ-1981 keyboard for an additional fee,
In the package:
Power cable,
2 GB USB flash drive,
Dimensions: 29.5 x 42.5 x 14.5 cm,
Weight: 11 KG,
Has current inspection and is ready for operation,
Technical Passport (Service Report) issued and valid for 12 months,
Warranty:
6 months on domestic market (Poland),
3 months on international market,
Possibility of extending the warranty for an additional fee up to 12 months for domestic and international markets,
Warranty period does not cover bulbs and other light sources installed in the device,
Financing options (Poland only): Installments, Leasing, Loan,

