



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Patient warmer MED WARM W-300 (Reconditioned)**

Refurbished (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
Manufactured in Turkey, (2015),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 350 W,  
System for patient warming  
MEDWARM  
is designed to prevent patient body cooling (hypothermia) by stabilizing or raising peripheral temperature. The modern system is intended for use with patients in operating rooms and post-anesthesia care units, as well as neonatal intensive care units, obstetric, orthopedic departments and ambulance services,  
Uses modern contact heating technology – heating material are carbon fibers,  
Safety:  
Precise, independently operating temperature sensors – each mat/blanket contains 10 temperature sensors,  
Dual control – each system is controlled by 2 independently operating microprocessors - one placed in the control unit, and the other in the mat/blanket,  
Independent safety shutoff - cuts power in case of overheating,  
System of audible and visual alarms,  
Low and safe voltage for heating elements: 24 V,  
Mattresses and blankets are radiolucent for RTG radiation – possibility to perform intraoperative radiographs,  
Patient comfort:  
Mattresses and blankets made of soft, flexible and lightweight material,  
Mattresses with memory foam layer reduce surface pressure on the patient's body,  
Temperature adjustment ranges:  
From 30°C to 40°C with 0.1°C steps,  
Built-in alarms:  
Temperature deviation / high temperature / disconnection of heating element / underheating (with mute option),  
Power supply loss alarm,  
Includes:  
Power cord,  
Heating mattress IM-120MS,  
Mattress dimensions: 120 x 50 x 3 cm,  
Total set weight: 10 KG,  
Has a current inspection and is ready for use,  
Technical Passport (Service Report) issued and valid for 12 months,  
Warranty:  
6 months on domestic market (Poland),  
3 months on international market,  
Option to extend warranty for an additional 12 months on both domestic and international markets,  
Financing option (Poland only): Installments, Leasing, Loan,

