



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SMITHS LEVEL 1 EQUATOR EQ-5000 Patient Warmer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Made in the U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Convection heating system
EQUATOR
provides controlled convection heating in a compact enclosure,
Temperature regulation electronically controlled using a thermoelectric sensor,
Can prevent all clinical effects of hypothermia, improve postoperative comfort of patients and shorten the length of stay in the post-anesthesia care unit,
The device can be free-standing or mounted on an IV stand, bed arm, or optional wheeled cart,
Effective temperatures:
High: 44°C +/- 1°C,
Medium: 40°C +/- 1°C,
Low: 36°C +/- 1°C,
Time to change the mean temperature from 20°C to 36°C: 7 minutes,
Equipment mounted on a dedicated mobile cart,
In the set:
Power cord,
New biological filter,
Air supply tube for warm air,
Dimensions: 24 x 30 x 26 cm,
Weight: 10 KG,
Has a current inspection and is ready for use,
issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Option to extend warranty for an additional 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

