



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET ALPHASTAR PLUS 1132.13A2 operating table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real product photos,
- German production,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in batteries 12 V / 17 Ah,
- Frequency: 50/60 Hz,
- Power: 400 VA,
- New upholstery (memory visco) with medical certificate in black,
- The entire 4-segment tabletop is radiolucent. The table features a REVERS function (inverted patient position),
- Mobile universal operating table with electrohydraulic drive. Emergency control panel on the table column,
- Adjustable in many planes using a wired remote control or via the panel on the column:
 - ON/OFF,
 - LOCK/UNLOCK - locking and unlocking the table,
 - NORMAL/REV - patient orientation on the operating table,
 - TREND/REV - changing the patient's lying direction relative to the table column: Trendelenburg and anti-Trendelenburg position within the range of $\pm 30^\circ$,
 - The table is equipped with an auxiliary drive controlled by a handheld remote control, allowing forward and reverse travel,
 - UP/DOWN - height of the operating surface from 68.5 to 107 cm (up/down),
 - TILT L/TILT R - lateral tilt in the range of 20° ,
 - BACK UP/BACK DN - adjustment of the back section position,
 - LEG UP/LEG DN - moving the connected leg segment: up/down,
 - FLEX/REFLEX - bending the tabletop in half in both directions: roof position / V position (Beach chair),
 - LEVEL - reset: return to the horizontal position of the bed,
- Manually adjustable:
 - 2-segment leg rest,
 - Headrest tilt,
- The base housing is made of durable plastic,
- The table has X-RAY transparency,
- Equipped with rails for mounting additional equipment,
- Table adapted to C-arm,
- Mobile, on double wheels
- Manually removable headrest and foot section,
- Included in the set:
 - Power cord,
 - MAQUET wired remote control,
- The table can be equipped with other accessories from our offer for an additional fee, examples are shown in the last photos,
- Dimensions:
 - Headrest: 52 x 30 cm,
 - Back support: 52 x 30 cm,
 - Bed segment: 52 x 98 cm,
 - Leg rest segment: 23 x 75 cm,
- Maximum total load capacity: 450 KG,
- Weight: 325 KG,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months to the international market,
 - Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,
- Financing options (Poland only): Installments, Leasing, Loan,

