



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Washer-disinfector PROHS WD 10 DIN D-10-SC/2 (Ex-Display)

New (on display),

Actual product photos,

Manufacturing: Portuguese, (2017),

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 400 W,

PROHS WD 10

with manually operated double doors for surgical departments and smaller sterilization departments with six washing levels and a connection for washing and disinfecting cylindrical instruments,

Features:

WIDE BASKET SELECTION

- a wide range of baskets allows washing and disinfection as well as drying even the most demanding instruments, including cylindrical shapes,

FAST CYCLE (SC) 50% SHORTER CYCLE

- additional tanks allow preheating water to speed up the washing phase. Total washing time is reduced by up to 50%,

COMPLIANT WITH ISO 15883 1/2

- The PROHS washer-disinfector was designed in accordance with ISO EN-ISO 15883-1/2,

FULLY AUTOMATED

- 20 standard programs. Some programs can be customized with individual parameters. The user has the option to choose a program suitable for their needs,

HEPA DRYING H14

- all instruments are thoroughly dried inside and out at the end of the cycle. The drying time and temperature are adjustable via the front control panel,

Technical specification:

Water pump capacity 626 L/min.,

Hot air drying with HEPA H14 filter - 150 m³/h,

Power consumption - 12 kW (13 kW SC model with fast cycle),

Touch color 3.5" control panel,

USB Port,

Two spray-wash arms (top and bottom of the chamber). Additional arms on each,

Built-in printer,

Vapor condenser,

Drawer for 4 chemistry canisters with a capacity of 5 L each,

Included:

4 pcs - dedicated footwear baskets,

3 pcs - dedicated tool baskets,

4 pcs - dedicated chemical canisters (capacity),

Vapor condenser,

Dimensions:

Chamber: 56 x 58 x 78 cm (capacity: 265 L),

Overall: 65 x 67 x 185 cm,

Weight: 280 KG,

Has current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months on domestic market (Poland),

3 months on international market,

Option to extend warranty up to 12 months for domestic and international markets at an additional cost,

Financing available (Poland only): Installments, Leasing, Loan,

