



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Linear USG probe MEDISON LH5-9ED (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

Made in Korea,

HL5-9ED probe

has frequency range from 5.00 MHz to 9.00 MHz,

Applications:

Thyroid,

Small parts of body,

Breasts,

Probe dedicated to ultrasound machines

:

Sonoace X4,

Sonoace 800,

Sonoace PICO,

Sonoace 8800,

Sonoace 6000,

Sonoace 5500,

Total weight: 1.5 kg,

Has current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend warranty for an additional 12 months for both domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

