



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



USG CURVED ARRAY HEAD MEDISON C3-7ED (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Korean production,
Transducer
Medison C3-7ED Curved Array
has a frequency range from 3.0 MHz to 7.0 MHz,
Applications:
Abdominal examination,
Fetal heart,
Obstetrics,
Gynecology,
Head/Probe dedicated to ultrasound machines
:
Medison SonoAce SA-8000SE,
Medison SonoAce SA-8000Live,
Medison SonoAce Pico,
Medison SonoAce X4,
Weight: 1.5 KG,
Has a current inspection and is ready for use,
Technical Passport issued (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend warranty for an additional fee up to 12 months for domestic and international market,
Financing available (Poland only): Installments, Leasing, Loan,

