



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



BERCHTOLD CHROMOPHARE D540 surgical lamp (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,

Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 240 VA,
Built-in bulbs: OSRMA 150 W / 22.8 V,
BERCHTOLD D540 surgical lamp

is intended to illuminate the operative field during surgical procedures. The lamp provides a high light intensity and uniform distribution on the operating field,

Specification:
Illuminance (at a distance of 1 m): 145 lux,
Range of adjustment of the illumination field: From 17 to 28 cm,
Color temperature of the light source (DIN 5035): 4500 K,
Electronic brightness regulation: From 65 to 145 lux,
Color Rendering Index (CRI): 94,
Height adjustment of the lamp head within 89 cm,
Free height adjustment, tilt of the head in the vertical and horizontal planes,
2-segment articulated arm with rotation adjustment range: 360°,
Built-in dual clamps for lamp positioning,
Autoclavable lamp holder,

Included in the set:
Ceiling cover for lamp mounting,
Minimum ceiling height for lamp mounting is 290 / 300 cm,
Dimensions:
Diameter of the body without railing: 49 cm,
Diameter of the body with railing: 54 cm,
Total length of the 2-segment arm: 189 cm,
Has a current inspection and is ready for operation,
Technical Passport (Service Report) issued and valid for 12 months,

Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend warranty with surcharge up to 12 months for domestic and international markets,
Warranty period does not cover bulbs and other light sources installed in the device,
Financing options (Poland only): Installments, Leasing, Loan,

