



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **SCHMITZ treatment lamp (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- LED bulb: 12V/50W,
- The lamp's exceptionally lightweight design allows for precise adjustment of the desired working position both vertically and horizontally. Instant ignition and light stabilization are achieved thanks to the use of an electronic converter.
- Elegant design and high-quality, durable materials. Its lightweight construction makes it suitable for use both for examinations in medical offices and as a user-friendly, individual light source in patient rooms. The lamp features a dedicated mount for various types of mounting rails with different diameters.
- Included in the set:
  - Attaching the lamp to the rail,
- Dimensions:
  - Mounting arm length: 78 cm,
- Weight: 2 KG,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Guarantee:
  - 6 months for the domestic market (Poland),
  - 3 months to the international market,
  - Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,
  - The warranty period does not cover light bulbs and other light sources installed in the device,
- Financing options (Poland only): Installments, Leasing, Loan,

