



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **NIDEK NT-2000 non-contact tonometer (Reconditioned)**

Refurbished (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
Made in Japan,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 90 VA,  
Software: 2.1,  
Technical specification:  
Measurement range: 1-60 mmHg,  
Built-in CRT 5" display,  
Working distance 11 mm (from air nozzle to cornea),  
Automatic pupil alignment,  
Built-in thermal printer (paper width 56 mm),  
Electrical eyelid lift,  
Automatic safety lock to prevent contact with the patient eye,  
Computer interface for data transmission,  
Manual chin regulation within 7 cm,  
Menu in English,  
In the set:  
User manual in English (PDF),  
Power cable,  
Dimensions:  
Overall: 26 x 52 x 45 cm,  
Weight: 19 kg,  
Category 02,  
Has a valid inspection and is ready for operation,  
Technical Passport issued (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Warranty extension available for an additional charge up to 12 months for both domestic and international markets,  
Financing available (Poland only): Installments, Leasing, Loan,

