



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER OXYLOG 3000 ambulance ventilator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Dutch production,
- Power supply: 230 V,
- Emergency power supply: via 1 new built-in battery 11.1 V / 7.8 Ah,
- Frequency: 50/60 Hz,
- **The Oxylog® 3000** is a time-variable, constant-flow, pressure-controlled respirator used as a rescue and transport respirator, providing volumes of 50 ml and above.
- **Scope of application:**
 - In transport in ambulances or helicopters,
 - In emergency departments, in recovery rooms,
 - When transferring patients requiring artificial respiration within a hospital or clinic,
- Built-in electroluminescent (EL) display,
- Displaying data in numerical and curved form,
- **Adjustable ventilator functions - setting values:**
 - Ventilation modes: IPPV, SIMV/ASB, CPAP/ASB, BIPAP/ASB,
 - Ventilation frequency (Freq) 2-60 1/min,
 - Tidal volume (Vt) 50 - 2000 mL,
 - Inhalation and exhalation time ratio 1:4 to 3:1,
 - Inspiratory time (Tinsp): 0.2 to 10 s,
 - Maximum airway pressure (Pmax) 20-60 mbar,
 - Inspiratory pressure (Pinsp): PEEP + 3 mbar to 55 mbar,
 - Mixing knob: 40 - 100% O₂,
 - Trigger sensitivity (FlowTrigger): 3 to 15 L/min,
 - Pressure booster ASB: 0 to 35 mbar
 - Pressure support rise time: slow, standard, fast,
- Feed gas: O₂,
- Power supply: from the central network or oxygen cylinder,
- Visual and audible security alarms,
- **Possibility to select menu in: GER, ENG, FR, ESP, PT, IT, NED, DUN, NOR, PL,**
- **Included in the set:**
 - User manual in Polish (PDF),
 - Power cord with power supply,
 - Patient hose system,
 - Gas supply hose with customer-selectable ends (AGA or DIN),
- Dimensions: 28.5 x 18 x 19 cm,
- Weight: 10 KG,
- Cat 02,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

