



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Ophthalmic instrument table SCHUMO SE-2 (Reconditioned)**

Refurbished (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

Polish production, (2015),

Power supply: 230 V,

Frequency: 50/60 Hz,

Table SCHUMO SE-2

is intended to cooperate with diagnostic ophthalmic devices,

Features a stable, durable and scratch-resistant base with adjustable aluminium feet, on which an electric column actuator from the Swiss company SCHUMO with a wide height adjustment range is mounted,

The actuator is mounted with a top measuring 78 x 43 cm,

Height adjustment range: 64 - 78 cm,

In the set:

Power cable,

Top dimensions: 78 x 43 cm,

Maximum load: 80 KG,

Weight: 40 KG,

Has a current inspection and is ready for operation,

Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Warranty extension available for an additional charge up to 12 months for both domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

