



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Aparat EKG NIHON KOHDEN ECG-2350 (Reconditioned)

Refurbished (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Manufactured in PRC, (2014),

Power supply: 230 V,

Backup power: via a new built-in battery,

Frequency: 50/60 Hz,

ECG machine

Nihon Kohden ECG-2350

is a 12-channel ECG machine with a 7-inch color LCD screen that simultaneously displays 12-lead ECG traces and data,

Internal memory can store the last 3 minutes of all ECG tracings. Large memory capacity ensures that essential information is not lost and also helps in identifying sudden arrhythmias,

Internal memory of the Nihon Kohden ECG device holds up to 400 ECG records, and on one SD memory card you can store about 3000 ECG files. ECG files (in DICOM or ECTP format) can also be transferred to a PC over a LAN network. The Nihon Kohden ECG machine allows viewing and storing data without paper,

Specification:

7-inch color LCD display enables easier simultaneous display of ECG traces and data from 12 leads,

Adjustable display improves visibility and provides smooth operation,

Ability to choose menu language: ENG, GER, FR, ESP, IT, RUS,

Mounted on a dedicated mobile cart,

In the box:

Power cable,

Dedicated cart for the machine,

KH-100D swiveling arm for cabling,

10-lead ECG cable,

Dimensions: 45 x 60 x 150 cm,

Total weight: 25 kg,

Has a current inspection and is ready for operation,

Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Warranty extension available for an additional fee up to 12 months for the domestic and international markets,

Financing options available (only Poland): Installments, Leasing, Loan,

