



PRAISTON Sp. z o.o.
ul. Górska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Ultrasonic disinfectant TROPHON EPR (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Manufactured in Australia, (2015),
Power supply: 230 V,
Frequency: 50/60 Hz,
Automatic high-level disinfection system for ultrasonic probes (HLD). Uses unique hydrogen peroxide mist technology (H₂O₂) that penetrates cracks and irregularities,
Fast (7-minute) process, ease of use, process quality assurance, environmentally friendly, no user exposure to harmful chemicals,
In the set:
Dedicated transport cart,
User manual in English, French (PDF),
Power cord,
Option to equip with consumables for an additional fee,
Dimensions:
Device: 34.5 x 34.5 x 49 cm,
Total: 50 x 60 x 140 cm,
Total weight: 38 KG,
Has current inspection and is ready for operation,
Technical Passport (Service Report) issued for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for international market,
Option to extend warranty for up to 12 months for domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

