



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Dr Mach 130F surgical-operating lamp (ceiling) (Reconditioned)

Refurbished (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

German production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Lamp: 50 W 24 V,

Dr Mach 130 F

is a halogen surgical operating lamp for wall mounting,

The head is equipped with a halogen lamp 22.8 V / 50 W,

The lamp produces a maximum light intensity of 50,000 lx,

Excellent color rendering,

Easy head position adjustment,

Ability to adjust light focus by rotating the holder,

Autoclavable holder,

Mounted on a jointed arm with total length about 155 cm,

Specification:

Color temperature: about 4300 K,

Color rendering index: 96 Ra,

Illumination field: Ø 14 - 25 cm,

Working distance: 70 - 140 cm,

Dimensions:

Head: 34 x 27 x 30 cm,

Total length of the 2-segment arm: 147 cm,

Weight: 14 KG,

Minimum ceiling height for lamp installation: 3 m (possibility to modify lamp for ceiling height of 2.5 m),

Includes ceiling cover plate,

Has current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months on the domestic market (Poland),

3 months on the international market,

Possibility to extend warranty for an extra charge to 12 months for domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

