



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Infusion fluid warmer LEVEL 1 SYSTEM 1000 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in U.K., (2019),
Power supply: 230 V,
Frequency: 50/60 Hz,
Smith Medical Level 1 H-1000
infusion fluid warmer that enables rapid, linear heating of IV fluids administered to the patient,
The fluid warmer provides rapid fluid flow to warm the patient during surgery. H-1000 uses disposable sets that allow blood and IV fluids to flow through the system and deliver fluids at a temperature of about 40°C,
Dimensions: 51 x 170 cm,
Weight: 10 KG,
Has current inspection and is ready for use,
Technical Passport issued (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend warranty for an additional 12 months for both domestic and international markets,
Financing options (Poland only): Instalments, Leasing, Loan,

