



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER FABIUS CE anesthesia machine (ventilator) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in batteries 12 V / 3.2 Ah,
- Frequency: 50/60 Hz,
- Software: 02.80,
- It has two ports for vaporizers with commonly used Dräger DW-2000 connector systems, **vaporizers not included - can be retrofitted for an additional fee,**
- Ventilator driven by an electric motor (electronically controlled),
- The airway monitor with control unit is a device for monitoring and measuring ventilation parameters during anesthesia,
- The monitor is equipped with a control unit that allows you to set and change all ventilation parameters listed below,
- **Gas supply:**
 - 2 x O₂,
 - 2 x N₂O,
 - AIR,
- Ventilator measurement functions:
 - Maximum pressure (P_{max}) 10 to 70 hPa,
 - Tidal volume (VT) 50 to 1400 mL,
 - Respiratory rate (f_{IPPV}) 6 to 60 L/min,
 - Inhale/Exhale Ratio,
 - (T_i : T_e),
 - 1:3 to 2:1,
 - Inspiratory pause (T_{ip} : T_i) 5 to 50%,
 - Positive End Expiratory Pressure (PEEP) 0 to 15 hPa,
 - Smooth manual gas flow regulation – O₂, AIR, N₂O, from 0.0 to 12L/min, (manometers),
- **Safety features:**
 - Monitoring P, V, FiO₂,
 - Automatic O₂ deficiency alarm and N₂O cut-off,
 - It has a built-in S-ORC (Sensitive Oxygen Ratio Controller), an oxygen ratio regulator with increased sensitivity (a control device ensuring a minimum O₂ concentration of 21 Vol.%),
 - Compact patient breathing system with oxygen sensor, placed on a movable arm,
 - Built-in electronic display with O₂, AIR, N₂O gas measurement function, as well as a control panel and ventilator function display,
 - Clear and easy to use interface,
 - Mobile device on 4 wheels, 2 of which with brakes,
- **Equipment:**
 - Extendable worktop with dimensions: 32 x 39 cm,
 - 2 drawers with dimensions: 32 x 32 x 12 cm,
 - Strips for mounting additional equipment,
 - 2 gas cylinder holders,
- **Possibility to select menu in: ENG, GER, FR, ESP, PT, IT,**
- Anesthesia machine equipped with a new O₂ sensor
- **Included in the set:**
 - DRAGER FABIUS CE operating instructions in Polish (PDF),
 - Gas supply hoses with ends to be chosen by the customer (AGA or DIN),
 - The patient's respiratory system,
 - Sampling line for gas measurement,
- **Possibility of equipping the device with a patient monitor/cardiometer from our offer - for an additional fee,**
- **Possibility of retrofitting vital signs monitoring cabling to the patient monitor/cardiometer - for an additional fee,**
- Dimensions: 70 x 102 x 110 cm
- Weight: 115 KG,
- Cat 05,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **The price of the device includes delivery, installation of the equipment and training,**
- **In the domestic market (Poland), if the device is installed within 30 days from the date of delivery of the equipment, the installation and training are free of charge, after exceeding this deadline these services will be charged,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

