



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Anesthesiologist's assistant METRO LIFELINE (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Mexican manufacture,
Mobile medical assistant
METRO LIFELINE
with extendable drawers can be mainly used in dental, surgical and procedural offices,
Technical specification:
Body, bumpers and handle made of plastic,
Acrylic covers on the top and front of the cart,
Includes tilting containers on the side panel,
IV pole with height adjustment,
Side mounting for cylinders,
Mobile top shelf for patient monitoring station,
Dimensions:
Patient monitoring station shelf: 46 x 35 cm,
2x retractable top drawer: 52 x 42 x 8 cm,
Retractable lower drawer: 52 x 42 x 16 cm,
Overall dimensions: 90 x 55 x 136 cm,
Weight: 50 KG,
Has current inspection and ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months on the domestic market (Poland),
3 months on the international market,
Option to extend warranty for a surcharge up to 12 months on both domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

