



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Karl Kaps KP 3000 Colposcope (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 150 VA,

Light source: 12 V / 100 W halogen lamp,
Colposcope application
KARL KAPS KP 3000

in clinical gynecology increases the accuracy of diagnoses and procedures performed. The KP 3000 colposcope was designed to help the doctor diagnose and treat,

High quality and comfort for physicians, the colposcope combines high-quality apochromatic optics and bright field illumination with ease and precision of maneuvering the colposcope, which results in high work ergonomics and great precision of procedures,

Specification:

Binocular 45° angle (f = 125 mm),

Objective: f = 300 mm,

Magnification adjustment: x0.63 / x1.0 / x1.6,

Eyepiece: WF 20 x V (with diopter adjustment),

Raised green filter,

Lighting: LEDone module, smoothly adjustable light intensity,

Manual height adjustment within 35 cm (head of microscope at heights 94 - 129 cm),

Built-in light source with free-adjustable intensity,

Overall dimensions: 50 x 50 x 94/129 cm,

Total weight: 55 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Option to extend warranty for a surcharge up to 12 months for domestic and international markets,

Warranty period does not include bulbs and other light sources installed in the device,

Financing options (Poland only): Installments, Leasing, Loan,

