



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Ophthalmoscope CARL ZEISS HO 110 (Reconditioned)**

Refurbished (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Germany,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Light source: 6 V / 5 W bulb,  
Optics and construction: High-quality optics and Build - sealed optical system,  
Equipped with a revolver full filter set:  
Blue,  
Green,  
Filament,  
Full regulation of correction using the revolver knob,  
Free adjustment of aperture shapes,  
Included:  
Power cord,  
Power supply,  
Total weight of the set 2 KG,;  
Has current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility to extend the warranty for an additional fee up to 12 months for both domestic and international markets,  
Warranty period does not cover bulbs and other light sources installed in the device,  
Financing available (Poland only): Installments, Leasing, Loan,

