



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Infusion pump (syringe) ALARIS GH (Reconditioned)**

Refurbished (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

Manufactured in U.K., (2017),

Power supply: 230 V,

Emergency power: via 1 new built-in Ni-Mh 7.2 V / 2.7 Ah battery,

Frequency: 50/60 Hz,

Power: 20 VA,

Alaris GH

is a multifunction syringe pump, used in intensive medical care as well as for other general medical applications,

Works with a wide range of standard disposable tip-lower-lock syringes. Compatible with syringes from 5 to 50 ml (e.g., B. Braun Omnifix, DB Plastipak, Monoject),

Features:

Large, readable graphical LCD display,

Three operating modes:

Injection rate,

Injection rate + infusion volume to be delivered,

Volume to be delivered + delivery time (automatic calculation of delivery speed),

KVO mode (keep-vein-open; maintaining catheter patency) with user-adjustable delivery rate in the range not less than 0.1 - 2.0 ml/h,

Drug library containing at least 100 drug names accessible from the pump keypad,

Automatic recognition of syringe size,

Delivery rate: from 0.1 to 1200 ml/h,

Memory for at least 1500 events,

RS232 connector,

Infrared communication port,

Convenient carrying handle,

Bolus delivery rate:

Syringe 5 ml - 10 - 150 ml/hour,

Syringe 10 ml - 10 - 300 ml/hour,

Syringe 20 ml - 10 - 600 ml/hour,

Syringe 30 ml - 10 - 900 ml/hour,

Syringe 50 ml - 10 - 1200 ml/hour,

Ability to choose menu language: ENG, NOR, GER, FR, DUT, SWE, ESP, PT, RUM, CS, PL, HUN, SVK, RUS, TUR, SERB, GREEK,

Included:

User manual in English (PDF),

Power cord,

Dimensions: 31 x 20 x 12.1 cm,

Weight: 2.7 KG,

Has current inspection and is ready for use,

Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for domestic market (Poland),

3 months for international market,

Option to extend warranty for an additional fee up to 12 months for domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

