



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Analytical scale SARTORIUS BL210S (Reconditioned)**

Refurbished (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

Manufacture: German,

Power supply: 230 V,

Frequency: 50/60 Hz,

The analytical balance is an extremely useful tool in laboratories, enabling precise and repeatable mass measurements of various substances,

Thanks to the advanced mechanism and modern electronics, analytical balances offer excellent accuracy and fast measurement speeds,

An important element of the analytical balance is the internal calibration system, which guarantees stability and precision of measurements,

even under intensive use. This ensures scientists and analysts can be confident that results obtained with the analytical balance are reliable

and meet expectations.

Specification:

Maximum load: 210 g,

Divisor (readability): 0.1 mg,

Linearity:  $<\pm 0.2$  mg,

Repeatability:  $<\pm 0.1$  mg,

Response time: 3 seconds (approx.),

Adaptation to ambient conditions: by selecting 1 of 4 optimized filtration levels,

External calibration mass: 200 g (class E2),

In the package:

Power supply,

Dimensions:

Draught: 8 cm,

Weighing chamber height: 20 cm,

Total: 18.9 x 25.1 x 29.9 cm

Weight: 3 KG,

Has current inspection and is ready to use,

Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Warranty extension available for an additional fee up to 12 months for both domestic and international markets

Financing options (only Poland): Installments, Leasing, Loan,

