



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MELAG Autoclave Type 17 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 1310 W,
Melag 17 is a steam sterilizer of class N,
Equipped with 2 operating programs:
Working pressure 1 bar - time 50 min - temperature 121° C,
Working pressure 2 bar - time 30 min - temperature 134 ° C,
On the lid of the autoclave there is a place to set down the trays,
Built-in manometer and thermometer,
Water discharge connection on the front of the device,
Pressure release button,
Indicators showing autoclave startup and heating,
All steam paths made of stainless materials,
Front feet with height adjustment,
Specification:
Chamber capacity: 11 L,
Water consumption during sterilization: 300 ml,
Maximum load: 3 kg of instruments including trays or 0.2 kg textiles,
Control panel in German,
Included in package:
Perforated tray with dimensions 40 x 14 x 2 cm (2 pcs.),
Tray handle,
Power cord,
Dimensions:
Overall: 63 x 43.5 x 35.5 cm
Chamber: Ø 17.5 x 48 cm,
Weight: 20 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend warranty for an additional charge up to 12 months for both domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan.



