



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE TRANSPORTER PRO patient monitor / cardiomonitor (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Made in USA,
Power supply: 230 V,
Backup power: via 2 new built-in batteries 11.1 V / 5 Ah (one battery included) - possibility to equip with a second one for an additional fee,
Backup power for the measurement module via a built-in new battery
Frequency: 50/60 Hz,
Cardi monitor / patient monitor
GE TRANSPORTER PRO
with built-in measurement module is a cardiomonitor monitoring multiple vital functions,
Ensuring continuous monitoring of vital signs during intra-hospital transport (ICU, operating room, ER, cardiology department),
Applications:
Adult patients,
Pediatrics,
Neonatal departments,
Specification:
Type: Active LCD matrix (TFT),
Diagonal: 10.4 inches (26.4 cm),
Resolution: 640 x 480 pixels (VGA),
Number of traces: Up to 6 leads simultaneously,
Roll speed: 6.25 / 12.5 / 25 mm/s,
Alarm system:
Visual: Alarm indicator integrated into the housing and parameter backlighting on screen,
Audible: Built-in speaker. Three warning levels (different tones),
Functions: Visual and audible indications of parameter exceedances, technical alarms (low battery, signal loss from module),
Includes:
Touch screen operation, and a navigation joystick allow direct access to many functions with simple commands,
Convenient handle and durable housing facilitate transport,
No separate hard drive and fan - the monitor is compact, quiet, stable and boots quickly,
Navigation dial facilitates entering information and navigating the screen,
Acoustic and visual alarms,
Built-in additional mounting,
Device equipped with GE measurement module with cabling:
ECG - 10-lead ECG cable,
SPO2 - new clip sensor,
NIBP - new cuff with air tube,
TEMP/CO,
IBP (invasive blood pressure) P1/P3 | P2/P4,
DEFIB/SYNC,
Menu language options: ENG, NED, FR, GER, IT, PT, ESP, DU, HUN, NOR, PL, RUS, SWE, CS,
Included:
English user manual (PDF),
Power cord,
Charger,
10-lead ECG cable,
New SPO2 clip sensor,
New NIBP cuff with air tube,
Dimensions: 30 x 33 x 15 cm,
Weight: 5 kg,
Has a current inspection and is ready for use,
Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for international market,
Option to extend warranty with an additional fee up to 12 months for domestic and international markets,
Financing option (Poland only): installments, leasing, loan,

