



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Digital flowmeter DRAGER DIGIFLOW O2 (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production, (2015),
- Powered by a built-in 3.6 V battery,
- **The DRAGER DigiFlow** device is designed based on modern sensor technology and a digital display,
- Flow measurement and settings are more precise than standard flow meters, and the regulators operate in smaller steps,
- **Characteristics:**
 - Accurate flow measurements,
 - Easy operation,
 - Clear reading of values,
- **Technical specifications:**
 - Parameter: standard (0-16 L/min),
 - Inlet pressure: 270-600 kPa (2.7-6 bar),
 - Flow range: 0-16 L/min,
 - Display range: 0.3-19.5 L/min.,
 - Flow restriction: approximately 20 L/min.,
 - Measurement accuracy: $\pm 10\%$ of the reading value,
 - Display accuracy: ± 1 digit,
- **Included in the set:**
 - Hose with AGA / DIN connection (to be selected when purchasing),
- **Flow meter prepared for humidifier,**
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

