



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FARUM A36 Oxygen Reducer (humidifier with doser) (Reconditioned)

Refurbished (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Polish production,

Oxygen reducer FARUM A36

with doser and humidifier with a reusable container for disinfection,

Oxygen reducer after connection to a oxygen cylinder, is intended for pressure reduction, setting the required flow, and humidifying oxygen drawn from the cylinder. The pressure in the cylinder is indicated by a manometer mounted on the reducer. Flow regulation allows delivering to the patient from 0 to 15 L/min of humidified oxygen. At the same time, it is possible to obtain a humidified oxygen-nitrogen mixture by turning the air-suction ring,

Specification:

Supply pressure: 2 - 15 MPa,

Reduced pressure: 0.4+0.1 MPa,

Flow regulation range: 0 - 15 L/min,

Gas humidification: approx. 30% RH,

Threaded nut $\frac{3}{4}$ "

Dimensions: 19.5 x 10 x 31 cm,

Weight: 1.8 KG,

Has current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Option to extend the warranty for an additional fee up to 12 months for both domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

