



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



### **DRAGER D300 surgical lamp (ceiling) (Reconditioned)**

Refurbished (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 80 VA,  
Lamp: 24 V 50 W,  
Technical specification:  
Single-focus,  
Halogen light,  
Light intensity (1m) - 35,000 Lux,  
Dome with adjustable holder in multiple planes,  
Mounting arm, 3-part allowing full control over the lamp,  
Includes switch placed on the lamp arm,  
Minimum ceiling height for lamp mounting: 240 cm,  
Dimensions:  
Total length on a 2-segment arm: 210 cm,  
Has current inspection and is ready for use,  
Technical Passport issued (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Option to extend warranty for an additional 12 months for domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan,

