



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDAP BORA P 2060 surgical suction unit (Reconditioned)

Refurbished (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

German production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 180 W,

MEDAP BORA P 206

0 is a small surgical suction device for aspirating secretions and fluids during procedures. It is mainly used in treatment rooms and hospital emergency departments (ED),

Modular overfill protection system

- the suction device features an innovative, modular overfill protection system to perfectly meet the specific requirements of applications both in operating rooms and in intensive care units and specialty departments,

Features:

Smooth adjustment of the maximum vacuum value,

Clear indicator of the current vacuum value,

Quiet and continuous operation,

Specification:

Operating mode: continuous operation,

Suction capacity (free flow): ~42 l/min,

Maximum vacuum: -90 kPa (-900 mbar),

Vacuum regulation range: 0 to -1000 mbar,

Noise level: 50.5 dB(A),

Included in the package:

Power cord,

SERRES jar with 2 L capacity and disposable insert,

Dimensions: 55 x 57 x 108 cm,

Weight: 21 kg,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend warranty for an additional charge up to 12 months for both domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

