



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ALVO MEDICAL ETIUDA 4.02 operating table (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Polish production,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in 12 V / 12 Ah batteries,
- Frequency: 50/60 Hz,
- Power: 350 VA,
- 6-segment,
- Base made of stainless steel,
- **Upholstery with medical certificate in black,**
- **Electric adjustment of settings using a wired remote control:**
  - Selecting head orientation,
  - Table "0" function,
  - Top-bottom range: 75-108 cm,
  - Vertical roll - Trendelenburg / Anti-Trendelenburg position, +35° / -35°,
  - Lateral tilt: +20° / -20°,
- **Manual adjustment of settings:**
  - Manually adjustable leg segment -100° / +20°, via gas springs, divided and 90° sideways,
  - Manually adjustable headrest +35° / -20° via gas springs,
  - Table top movement range: 30 cm,
  - Manually adjustable backrest +80° / -20° via gas springs,
  - Manual adjustment of the lumbar section in the range of: 11 cm (crank),
  - Manual brake,
- **It has:**
  - X-RAY translucency,
  - Rails for mounting additional equipment,
- **Dimensions of the bed:**
  - Width: 50 cm,
  - Length: 210 cm,
  - Thickness: 6 cm,
- **Mobile table with built-in central brake,**
- **Included in the set:**
  - Power cord,
  - New ALVO wired remote control,
  - Crank for adjusting the lumbar segment,
- **Possibility of equipping the table with additional accessories visible in the last photos,**
- Dimensions:
  - Base dimensions: 58 x 99 cm,
  - Overall dimensions: 58 x 210 x 75/155 cm,
- Maximum load capacity: 250 KG,
- Weight: 240 KG,
- Cat 02,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months to the international market,**
  - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

