



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Ultrasound Machine LOGIQ F6 (Reconditioned)

Refurbished (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Manufactured in P.R.C. (2016),

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 400 VA,

LOGIQ F6 ultrasound system

enables the acquisition of excellent image quality across a wide range of studies. It was developed to optimize workflow and increase diagnostic confidence, which in turn provides patients with higher quality medical care.

The presented system is an innovative solution that will facilitate patient diagnostics,

Unique imaging quality provided by the LOGIQ F6 system is invaluable when making a medical diagnosis. A clear image helps detect abnormalities and even small pathological changes,

Applications:

Abdominal cavity,

Obstetric and gynecologic,

Small organs,

Musculoskeletal system,

Vascular and peripheral vessels,

Urological,

Pediatric,

Intracranial,

Cardiologic,

Basic imaging modes:

B-Mode (2D) - basic imaging,

M-Mode (motion) - standard functional mode,

Color Flow Doppler (CFM) - color flow imaging,

Power Doppler / PDI - power Doppler,

PW Doppler (Pulsed-Wave) - spectral Doppler,

Ability to select language menu: ENG, CHIN, NOR, GER, IT, FR, ESP, RUS, SWE, BR, FI, NED,

The ultrasound device is equipped with a SONY UP-D711MD thermal printer,

Included:

Linear transducer GE L6-12-RS,

Overall dimensions:

Weight: 56 KG,

Has current inspection and is ready for use,

Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Option to extend warranty for an additional 12 months for domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

