



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ZOLL M SERIES biphasic defibrillator (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

Made in the U.S.A.,

Power supply: 10 V / 29 V,

Emergency power: via 1 built-in 10 V / 2.5 Ah battery,

Frequency: 50/60 Hz,

Power: 130 W,

Zoll M Series

are light and portable devices combining defibrillation, ECG, patient status monitoring with printing and data recording capabilities,

The defibrillator is intended for use during any type of resuscitation, and its light weight, small dimensions and compact construction make it particularly suitable for use during transport,

Specification:

Color LCD screen with a diagonal of 5.5 inches,

ECG module,

Printer for 90 mm wide paper,

Slot for external PCMCIA memory cards,

Menu language: English,

Control panel keyboard in English,

ZOLL M SERIES

is a portable defibrillator with a built-in handle and a dedicated t bag for wiring,

In the set:

Polish-language user manual (PDF),

Power cable for vehicle,

4-lead ECG cable,

AED cable,

Dimensions: 36 x 47 x 22 cm,

Weight: 10 KG,

Category D157,

Has a current inspection and is ready for operation,

Technical Passport (Service Report) issued, valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Option to extend the warranty for a surcharge up to 12 months for the domestic and international markets,

Financing options (Poland only): installments, leasing, loan,

