



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## CO2 incubator HERAEUS HERACELL 240 (Reconditioned)

Product temporarily unavailable  
Refurbished (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
German manufacture,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 640 W,  
HERACELL 240

is the ideal CO2 incubator for cell and tissue culture. The copper interior of this CO2 incubator is perfect for controlling contamination factors, as it is naturally fungicidal and bactericidal, and also an excellent heat conductor. The Heracell 240 also features a built-in ContraCon 90°C disinfection procedure that effectively removes all contaminating organisms in a sterile environment,

Equipped with a microprocessor-based control and an easy-to-read digital interface that allows users to easily adjust settings and alarms, Tempered glass doors allow users to review the progress of their samples without disturbing chamber operation. The incubator also regenerates 5 times faster than other models when opening is necessary, in only 10 minutes at 37°C,

Humidity is more reliable using a tilted reservoir instead of a traditional basin, which eliminates the risk of transfer and contamination and includes a low-level warning sensor,

HERACELL 240 operates in a temperature range from +5°C above ambient to 55°C (+/- 0.5°C at 37°C). CO2 content in the chamber can be adjusted in the range 0 - 20% (+/- 0.1%). Relative humidity is maintained at 90% or 95% as set by the user,

Chamber capacity: 240 L,  
Keyboard on the control panel in English,  
Included:

English-language user manual (PDF),  
2 x tray,  
Power cord,

Dimensions:  
General: 82 x 78 x 94 cm,  
Chamber: 60 x 66 x 53 cm,  
Tray: 56 x 50 cm,  
Weight: 107 kg,

Has current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:

6 months for the domestic market (Poland),  
3 months for the international market,

Possibility to extend warranty up to 12 months for both domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan,

