



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MEDELA DOMINANT 50 surgical suction pump (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Swiss production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 150 W,
- **Dominant 50** is a high-quality surgical suction unit that guarantees maximum suction efficiency and quiet, reliable operation,
- **Product features:**
  - Maintenance-free low-speed piston pump,
  - Use in procedures requiring rapid suction of large amounts of fluids,
  - Recommended mainly for operating rooms, neurosurgery, plastic and general surgery,
  - Piston-cylinder technology ensures achieving the desired vacuum in a few seconds and high suction efficiency,
  - No vibrations thanks to low rotational speed,
  - Designed for continuous operation (24 hours a day),
  - Long life of the device,
  - Smooth suction power regulation,
  - Built-in pressure gauge,
  - The diaphragm vacuum regulator allows for precise vacuum setting,
  - A safety kit located on the side of the pump prevents overflow and reduces the risk of contamination,
  - The base has a rail for mounting jars and accessories,
  - Rail mounting using bones,
  - Possibility of using jars of different capacities,
- **Specification:**
  - Capacity: 50 L / min,
  - Smooth suction power regulation: max. 90 kPa (675 mm Hg),
  - Time to obtain vacuum from 0 to 85 kPa on a 2 L tank: 14 s,
- **Included in the set:**
  - User manual in: EN, DE, FR, IT, NL, SV, FI, NO, ES, PT (PDF),
  - 2 x autoclavable jar (121° C) SERRES, capacity 2 l, with disposable insert,
  - Suction hose 150 cm long,
  - Power cord,
- Dimensions:
  - Suction: 39 x 25 x 33 cm,
  - Overall: 48 x 48 x 108 cm,
- Weight: 24 kg,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months to the international market,**
  - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

