



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HERAEUS UT12 Laboratory Warmer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1.82 KW,
- **The HERAEUS T12 incubator** utilizes natural convection, making it suitable for thermal process applications such as heated storage, heat treatment, thermal testing, and convective drying processes at temperatures up to 250°C. The incubator is equipped with a precise microprocessor controller with a large, easy-to-use display.
- **Specification:**
 - Temperature range adjustment: Ambient temperature +10°C to 250°C,
 - Internal capacity: 110 liters,
 - Maximum temperature: 250°C,
 - Precise microprocessor controller,
 - Natural convection,
 - User-friendly operation,
 - Easy to clean and maintain - with rounded corners and removable shelf supports,
 - Corrosion-resistant internal construction made of stainless steel,
- **Temperature deviations:**
 - Spatial:
 - 70°C (±3.0°),
 - 150°C (±5.0°),
 - 250°C (±6.0°),
 - Temporary:
 - 70°C (± 0.3°),
 - 150°C (± 0.5°),
 - 250°C (± 0.7°),
- The incubator is equipped with one tray measuring 36 x 47 cm,
- Dimensions:
 - Total: 69.6 x 54 x 85 cm,
 - Internal: 49.6 x 38.6 x 57.5 cm,
- Weight: 55 KG,
- Cat 02,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

