



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Urology chair STOLTER FG-04 (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Polish production,
Power supply: 230 V,
Frequency: 50/60 Hz,
New upholstery with medical certificate in gray color,
Urology chair FG-04
- intended for conducting gynecological and urological examinations and procedures in healthcare facilities,
Construction made of steel profiles coated with powder coating,
Upholstery with medical certificate in gray,
Specification:
Electrical adjustment using foot controller:
Backrest tilt range: 66°,
Chair tilt: -12° / 0° / 20°,
Possible chair position settings:
Trendelenburg: 24°,
Anti-Trendelenburg: 12°,
Equipped with:
Stainless steel rails allowing attachment of additional equipment,
2x removable leg supports,
Removable stainless steel sink with dimensions: 32 x 26 x 10 cm,
Paper holder with width 33 cm,
Base of chair with a folding platform,
Additional sliding leg support segment,
Features:
Leg supports with wide-range adjustable,
Possibility to change the patient holder position,
Dimensions:
Backrest: 62 x 85 x 6 cm,
Seat: 62 x 37 x 6 cm,
Base: 60 x 100 cm,
Leg support segment: 60 x 45 cm,
Overall: 146 x 80 x 158 cm
Maximum load: 175 KG,
Weight: 75.5 KG,
Has current inspection and is ready for use,
Technical Passport (Service Report) issued and valid for 12 months,
Warranty:
6 months on the domestic market (Poland),
3 months on the international market,
Option to extend warranty for an additional charge up to 12 months for domestic and international markets,
Financing options (only Poland): Installments, Leasing, Loan,

