



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **PHILIPS EVERFLO OPI Oxygen Concentrator (Ex-Display)**

New (display),

Product actual photos,

Made in Mexico, (2025),

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 300 W,

Philips Respironics Everflo

is a lightweight and easy-to-use oxygen concentrator,

The device generates concentrated oxygen from ambient air drawn from the surroundings to deliver to a patient requiring low-flow oxygen therapy. Oxygen from the air is concentrated using a molecular sieve in a pressure swing adsorption process,

Quiet and economical,

Specification:

Equipped with a flow meter,

Smooth flow adjustment in the range of 0.5 - 5 L/min,

Oxygen concentration (at 5 LPM): 93%,

Noise level: < 40 dB,

Mobile on wheels, all with braking,

Included:

Polish language user manual (PDF),

Power cord,

6x new oxygen mask size L / 210 cm with drain,

Humidifier (bottle),

Dimensions: 38 x 25 x 59 cm,

Weight: 15 kg,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend the warranty for an additional charge up to 12 months for the domestic and international markets,

Financing available (Poland only): Installments, Leasing, Loan,

