



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RUSCH LIFE-AIR 1000 patient warmer (Reconditioned)

Refurbished (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

German production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 120 VA,

Hypothermia therapy system

LIFE - AIR 1000

was developed to reduce the effects of hypothermia, which is a common issue during and after a procedure,

Manual adjustment of the air temperature flowing around the patient using the control switch heater in the range to choose from:

27° C,

32° C +/- 3° C,

36° C +/- 3° C,

43° C +/- 3° C,

Low weight facilitates transport and preparation for work,

Mobile on double wheels,

In the kit:

Operating manual in English (PDF),

Power cable,

Dimensions:

Heater: 29 x 30.5 x 38/58 cm,

Overall: 53 x 53 x 75 cm,

Weight: 15.5 KG,

Has a current inspection and is ready for use,

Technical Passport (Service Report) issued for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility of extending the warranty for an additional charge to 12 months for the domestic and international market,

Financing options (Poland only): Installments, Leasing, Loan,



