



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Block heater YWR DIGITAL HEATBLOCK (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in USA,

Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 210 W,

The YWR DIGITAL HEATBLOCK heater

is designed for applications requiring repeatable results and excellent temperature stability. These multifunction devices are ideal for incubation, enzymatic reactions, immunological tests, melting and boiling point determination, and a wide range of other laboratory procedures,

Microprocessor PID control provides stability and temperature uniformity at  $\pm 0.2$  °C,

LED display shows current temperature and set temperature. Controllers allow setting the temperature in steps of  $\pm 0.1$  °C.,

Specification:

Temperature range: ambient temperature + 5°C to 120°C,

Stability:  $\pm 0.2$ °C at 37°C,

Stability:  $\pm 0.4$ °C at 60°C,

Achieving 100°C is possible in the range from 45 to 65 minutes (depends on block capacity),

Included in the set:

Power cord,

Dimensions: 20.3 x 46.7 x 8.9 cm,

Weight: 3.6 KG,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months in the domestic market (Poland),

3 months on the international market,

Warranty extension available for an additional charge to 12 months on the domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

